

# Northern Lights user manual




## Basic user guide – v6 Jan 2026

<a href="#"><u>Start-up</u></a>	<a href="#"><u>2</u></a>
<a href="#"><u>Daily QC</u></a>	<a href="#"><u>3</u></a>
<a href="#"><u>Setting up a new experiment</u></a>	<a href="#"><u>4</u></a>
<a href="#"><u>Sample Acquisition</u></a>	<a href="#"><u>8</u></a>
<a href="#"><u>Unmixing</u></a>	<a href="#"><u>10</u></a>
<a href="#"><u>Exporting data</u></a>	<a href="#"><u>11</u></a>
<a href="#"><u>Analysis</u></a>	<a href="#"><u>12</u></a>
<a href="#"><u>Standby wash</u></a>	<a href="#"><u>14</u></a>
<a href="#"><u>Shutdown</u></a>	<a href="#"><u>14</u></a>
<a href="#"><u>Waste tank full</u></a>	<a href="#"><u>15</u></a>
<a href="#"><u>Sheath tank empty</u></a>	<a href="#"><u>16</u></a>

*Click titles to return to this page*

## Start-up

NL has a 30 min warm up time. Timer starts after opening a template



Spectroflo icon

- 1) Check sheath / waste tank levels
- 2) Power on the cytometer. Power button on the left-hand side of the cytometer.

Windows credentials: NLuser; P: Welcome#1

- 3) Open the SpectroFlo software
- 4) Log in

Start typing your name in the username field; mouse-click to select user account  
Password: first letter of first name; entire last name; all lowercase


QC/startup account U: QC\_startup | P: CytekNL

### *Run water during warm up to flush fluidics*


- 1) Click the [acquisition] icon on the welcome screen
- 2) Select the [Default] experiment layout
- 3) Check cytometer is connected (green tick, bottom right)
- 4) Remove tube from the sample probe
- 5) Load a tube with 2.2mL milliQ water
- 6) Set the flow rate [high]
- 7) Click on [start ►] and run water as the cytometer warms up for at least 5 minutes.




when changing sample tubes, wait until the SIT flush is complete.



New default experiment



Welcome screen




Flow rate

## Daily QC


- 1) SpectroFlo QC Beads are in the fridge under the centrifuge
  - a. Check if there is already a 5 mL tube made up (need >80  $\mu$ L for a QC run)
  - b. *If need to prepare a new bead tube:*  
Add 300  $\mu$ L milliQ to a new 5mL glass tube  
Vortex bead bottle  
Add one drop to glass tube. Vortex glass tube.



- 2) Click on the QC tab in the blue ribbon, top right




- 3) Load 5 mL bead/milliQ tube on the cytometer
- 4) Check bead lot
  - a. Bead lot on the side of bead bottle
- 5) Click start
- 6) QC will take 3-5 minutes
  - a. After QC finishes, return bead bottle and bead/milliQ tube (if >100  $\mu$ L) to fridge




## New experiment

- 1) Click on the [acquisition] tab on the blue ribbon.



- 2) Close the start-up experiment.

(X) to right of experiment name; to left of "WORKSHEETS"



- 3) Select [New] experiment layout




- 4) Enter exp name (YYYYMMDD\_exp name)

- 5) Choose parameters


Expand library subgroup and double-click on fluorophore needed  
If not there, use the closest fluor



Click > to expand



- 6) Click [next] (bottom right)
- 7) Add a reference group (single stain controls)



## 8) Setup single stain controls

Set cells or beads.

Set [negative control] for fluors if these controls do not have an "internal negative population (ie, CD45 or cells that are 100% GFP+).

Define Unstained Control(s) for Autofluorescence Extraction

Name Unstained


Control Type Cells

Select if used beads AND cells for controls


☐ Define Additional Negative Control(s) for Spillover Calculation

Fluorescent Tag	Control Type	Label	Lot	Negative Control
BFP	Cells			
BB515	Cells			
APC	Cells			

## 9) Add a new group for the FMOs and full stain samples



## 10) Set names for the tubes



## 11) Click [next]

## 12) Set marker names

Click on the box at the experiment level (see below) to set the name for all tubes (in the image below, did not press [enter] yet, so it hasn't applied label to all.

Can also set labels at the group level

Will not apply to reference controls

Fluorescent Tags Groups Markers Keywords Acquisition				
Edit Lot				
Groups				
Name	BFP	BB515	APC	
20231115_testy mctestFace		CD45		
Reference Group				
Unstained (Cells)	N/A	N/A	N/A	
BFP (Cells)		N/A	N/A	
BB515 (Cells)	N/A		N/A	
APC (Cells)	N/A	N/A		
samples				
patient 001				

Experiment level

Enter label here to apply to all tubes


Group level

- 13) Click next
- 14) Keywords (can leave blank)
- 15) Click [next]
- 16) Set stopping event number, volume, or time.  
Again, can set at a global level, group level, or individual tubes

☐ Tube Specific User Setting    Experiment User Setting: CytekAssaySetting (Cytek) ▼

Name	Worksheet	Stopping Gate	Storage Gate	Events To Record
▼ 20231115_testy mctestFace	Default Raw Worksheet (Raw) ▼	All Events ▼	All Events ▼	5,000 ▼
▼ Reference Group	Default Raw Worksheet (Raw) ▼	All Events ▼	All Events ▼	5,000 ▼
Unstained (Cells)	Default Raw Worksheet (Raw) ▼	All Events ▼	All Events ▼	5,000 ▼
BFP (Cells)	Default Raw Worksheet (Raw) ▼	All Events ▼	All Events ▼	5,000 ▼
BB515 (Cells)	Default Raw Worksheet (Raw) ▼	All Events ▼	All Events ▼	5,000 ▼
APC (Cells)	Default Raw Worksheet (Raw) ▼	All Events ▼	All Events ▼	5,000 ▼
▼ samples	Default Raw Worksheet (Raw) ▼	All Events ▼	All Events ▼	5,000 ▼
patient 001	Default Raw Worksheet (Raw) ▼	All Events ▼	All Events ▼	5,000 ▼


- 17) Click [save and open]





## Sample acquisition

Final protocol checks


- Flow rate (can not change mid run; requires [pause])
- User settings (loads fluorescence and scatter gains and threshold)
-  DO NOT change individual channel gains. If signal too high or low, change via the “all channel %” setting.

Ballpark user-settings (derived from cytek assay settings)


	leuk (human)	Cell lines
FSC	50	20
SSC (violet / blue)	145 / 160	60 / 60
Threshold	FSC, 350,000	FSC, 350,000
V1	45	20
B1	1090	50
R1	80	20

Use the all channel % option to set V1, B1 and R1 to the numbers above. This will also set the other channels V2-V16, B2-B14, R2-R8 to the appropriate gain (relative to % change).

Changing ‘all channel %’, increases or decreases all the gains of that laser line by the number (percentage) entered. The % value is relative to the loaded user settings. So, if you save these settings, the % value will be set to 0, with any changes from then being applied to the newly saved gains.



Choosing predefined user settings



Setting FSC and SSC

GAIN THRESHOLD SIGNAL LASERS

FSC SSC SSC-B

50 145 161

Violet Blue Red

V3 V4 V5

46 69 55 40 55

V6 V7 V8 V9 V10

46 61 85 64 73

V11 V12 V13 V14 V15

63 57 51 61 83

V16

70

All Channels %: 0

Change all gains for the selected laser

GAIN THRESHOLD SIGNAL LASERS

Threshold Operator: ☒ Or ☐ And

Channel	Threshold
FSC	350,000


Threshold settings

- 1) Load unstained sample, this will be the first tube in the group.

*If you want, you can load a full stain tube to check if the fluorescence signal is on scale. Then switch back to the unstained to record data for that tube.*

- 2) Press [Start ▶]

Make sure that scatter and snake looks ok (ballpark: set lowest channel (of each laser near 0))



Scatter and snake for unstained PBMCs

- 3) Press [Record] to save data for that tube

Once finished, the selected tube will go down the line of the group. You can manually change tubes by left-clicking on them.




when changing sample tubes, wait until the SIT flush is complete.


- 4) Run/record all reference controls and FMO/full-stains.

## Unmixing

- 1) Click unmix



- 2) Set the unstained reference tube (beads or cells).  
Leave blank if reference tube has internal negative population.
- 3) Select AF as a fluorescence tag (bottom left).



- 4) Click [next]



- Set the scatter gate appropriately and then set the positive and negative populations for each fluorophore.

FSC/SSC: Tight selection around cells of interest.


Positive peak: Mode and greater, or brightest 500-1000 events

Negative peak: Mode and lower

You can manually change the dominant channel used for defining pos / neg by dragging the “error bar” on the snake graph to another channel.





- Select the [live unmix] button to apply unmixing to this experiment (the non reference group samples).



## Exporting data

Export files to open in FCS express or flowJo.

- Close experiment
  - Click [my experiments]
  - Right click experiment of interest
  - Click export
  - Save on D:\ in your folder
  - Raw and unmixed files exported.
- Open the unmixed file in analysis software




## Analysis

You can export the unmixed data to FCS express or plot data in spectroFlo.

- 1) Create a new unmixed worksheet

Clicking the (+) to the left of the [REPORTS]

Select [new unmixed]





- 2) Create plots and gates

To change the parameter, right click the plot axes

To set a gate to a plot, right click the plot > properties


The properties window also allows bi-ex limits, dot size.






## **Standby wash**

- 1) Prepare 2 tubes
  - a. 2.2 mL sodium hyperchlorite
  - b. 2.2 mL water
- 2) Click the third icon on the left (cytometer)
- 3) Click [clean flow cell]
- 4) Follow on-screen instructions  
/!\ except for reagent volumes and concentrations. These are incorrect.
- 5) Log out of SpectroFlo



## **Shutdown**

- 1) Click the third icon on the left (cytometer) > fluidics shutdown
- 2) Prepare 3 tubes (2.2 mL):
  - a. Bleach (1:2 bleach : water)
  - b. Water
  - c. Contrad (1:9 contrad : water)
  - d. Water, can refill the second tube.
- 3) Follow on screen instructions.  
/!\ except for reagent volumes and concentrations. These are incorrect.
- 4) When finished, software will instruct to power down cytometer (power button on side of machine)
- 5) Log out of SpectroFlo
- 6) Leave PC on, it will go into standby.





## Sheath and waste


### *Waste tank full*

Decontaminate the full waste tank

- 1) Unplug fluidics line on the waste tank lid (quick connect)
- 2) Unplug electronics cable from the cytometer body (black wire)
- 3) Add 350 mL sodium hypochlorite to decontaminate waste  
Bleach bottles pre-filled with 350mL bleach on trolley.  
/!\ Do NOT throw out hypochlorite bottles; they will be reused.  
Tilt the bleach bottle on a 45-degree angle to denote it being empty
- 4) Write date bleach added on lid (whiteboard marker)
- 5) Place waste tank on trolley.

*Connect empty waste tank to cytometer*

- 6) Put empty waste tank in the holder
- 7) Attach the fluidics and electronics lines.
- 8) Resume your cytometer run.



Waste tank in biocan holder  
on flow lab trolley





### Sheath tank empty

- 1) Push in metal tab on sheath tank to disconnect fluidics line from the sheath tank lid.
- 2) Remove empty sheath tank from the holder
- 3) install the full sheath tank in the holder and connect fluidics line.
- 4) Place empty sheath tank on trolley

Or:

- 5) Refill empty sheath tank. Use Satorius system in adjacent room.
- 6) Unscrew lid and fill tank with 3.8 L of milliQ water (2 x 1.9 L).  
Biocan holder available to aid in transporting between the chem prep room and 9-50.
- 7) Place the newly filled sheath tank on black shelves.

