

# Labelling Human Cells for Flow Cytometry: Direct Protocol

Please note – this protocol has been optimised for certain conditions and is meant as a guide only, you may need to alter some parameters to suit your own experiment

1. Wash cells at least once and resuspend to approx  $4 \times 10^6$  viable cells/ml in staining buffer.
2. Add purified human IgG to block Fc receptors. Use 5  $\mu$ l of 10 mg/ml stock per million cells. Incubate at room temperature for 30 minutes
3. Aliquot 50  $\mu$ l of Fc-blocked cell suspension into FACS tubes (Falcon 2008)
4. Add 10  $\mu$ l of each primary antibody and mix. For all primary antibodies the optimal concentration should be determined by titration. Each antibody must be directly labelled with a fluorochrome or biotin, and all antibodies that are to be used together in one tube must be labelled with *different* conjugates.
5. Incubate at room temperature for 30 minutes (*during incubation, set centrifuge to cool to 4°C and dilute streptavidin conjugate if required... see Step 7*)
6. Add 3 ml of ice-cold staining buffer, spin down, tip off supernatant (*allow to drain briefly after tipping off, and be consistent between tubes, as different amounts of liquid left in bottom of tube may affect subsequent steps*)
7. If biotinylated antibodies are used, then follow these steps. If all antibodies are directly labelled with fluorochromes, go straight to Step 13.
8. Dilute streptavidin-fluorochrome conjugate to recommended concentration in staining buffer
9. If using a cyanine-containing dye such as PEcy5, add human IgG to 1% and allow to pre-adsorb for at least 20 minutes at room temperature
10. Add 50  $\mu$ l of streptavidin conjugate to all tubes containing biotinylated mAb. Vortex for 3 seconds
11. Incubate on ice 30 minutes
12. Add 3 ml of ice-cold PBS/azide (only to those tubes that received the conjugate), spin down, tip off supernatant, vortex briefly.
13. Add 3 ml ice-cold PBS/azide to ALL tubes. Spin down, tip off supernatant
14. Resuspend in 200 – 500  $\mu$ l 1% paraformaldehyde (depending on number of cells) to fix
15. Store at 4°C in the dark until analysis

## Suggested controls\*:

1. Cells only
2. Isotype-matched negative control antibodies for each fluorochrome, or FMO controls
3. Compensation controls (ie, a brightly-staining antibody by itself) for each fluorochrome

**\*NOTE:** See Recommended Controls document for comprehensive details