

Flow Cytometry (murine): Polyclonal Antibody

Please note – this protocol has been optimised for certain conditions and is meant as a guide only, you may need to alter some parameters to suit your own experiment

1. Wash cells in PBS and resuspend to approx 4×10^6 viable cells/ml in PBS/BSA/azide
2. Add 5 μ l of murine gamma-globulin (at 10 mg/ml) per 10^6 cells, to block Fc receptors.
3. Incubate at room temperature for at least 20 minutes (during incubation, prepare secondary antibody - see Step 8). There is no need to wash out the gamma-globulin.
4. Aliquot 50 μ l of Fc-blocked cells into FACS tube (Falcon 2008)
5. Dilute primary antibody (purified IgG) to an appropriate concentration in PBS/BSA/azide. This will need to be determined by titration.
6. Add 10 μ l of primary antibody to cells, pipette up and down to mix, and incubate at room temperature for 30 minutes (during incubation, set centrifuge to cool to 4°C)
7. Wash with 3 ml ice-cold PBS/BSA/azide
8. Add 50 μ l of secondary antibody (eg anti-IgG (or anti-IgM) conjugated to a fluorochrome or biotin) to cells. *Secondary antibody should be diluted to an appropriate concentration in PBS/BSA/azide To block non-specific binding, add 2% normal mouse serum and 100 μ g/ml mouse gamma-globulin and pre-adsorb for 30 minutes prior to use.*
9. Incubate on ice for 30 minutes.
10. Wash with 3 ml ice-cold PBS/BSA/azide
11. Dilute streptavidin-conjugate to appropriate concentration in PBS/BSA/azide and add 50 μ l to cells. Incubate on ice for 30 minutes
12. Wash with 3 ml ice-cold PBS/BSA/azide then with 3 ml ice-cold PBS/azide (no BSA)
13. After final spin, tip off the supernatant, vortex, add 200 – 500 μ l of 1% paraformaldehyde and vortex briefly again
14. Store in dark at 4°C until analysis (up to 2 weeks)
15. Note: PBS/BSA/azide = PBS containing 1% BSA and 0.04% sodium azide

Suggested controls*:

1. Cells only
2. Isotype-matched negative control antibodies for each fluorochrome, or FMO controls
3. Compensation controls (ie, a brightly-staining antibody by itself) for each fluorochrome

***NOTE:** See Recommended Controls document for comprehensive details