

Recommended controls

You need to include adequate controls for your experiment to have confidence in drawing conclusions from, and interpreting the data.




Consequences:

Conclusions or hypothesis founded on false positive/negative data.

The data acquired cannot be used for publication – so, you need to repeat the experiment.













Single colour experiments

Controls for

Unstained cells		Setting voltage of your fluorophore (auto-fluorescence) Shows where negative cells reside in your data Setting FSC and SSC of your cell type
Isotype control (antibody-based flow)		Non-specific binding of your fluorescent-antibody
Experimental control (e.g. ± compound)		Cell population FSC and SSC changes caused by experiment Changes in auto-fluorescence

Multi-colour experiments

All of the controls mentioned above, plus:

Compensation controls	<table border="0" style="display: inline-table; margin-right: 10px;"> <tr><td style="text-align: center;">FITC</td></tr> <tr><td style="text-align: center;"></td></tr> </table> <table border="0" style="display: inline-table; margin-right: 10px;"> <tr><td style="text-align: center;">PE</td></tr> <tr><td style="text-align: center;"></td></tr> </table> <table border="0" style="display: inline-table; margin-right: 10px;"> <tr><td style="text-align: center;">APC</td></tr> <tr><td style="text-align: center;"></td></tr> </table> <table border="0" style="display: inline-table;"> <tr><td style="text-align: center;">Cy7</td></tr> <tr><td style="text-align: center;"></td></tr> </table>	FITC		PE		APC		Cy7		Single-stained tubes for every fluorophore in your assay Needs to be performed EVERY; yes, EVERY flow acquisition
FITC										
										
PE										
										
APC										
										
Cy7										
										

Controls for

Spectral overlap of your fluorophores into channels other than its own
Visualise spectral spread after applying compensation

FMO controls	<table border="1" style="border-collapse: collapse;"> <thead> <tr> <th></th> <th style="background-color: #d9ead3;">FITC</th> <th style="background-color: #f4cccc;">PE</th> <th style="background-color: #f4cccc;">APC</th> <th style="background-color: #f4cccc;">Cy7</th> </tr> </thead> <tbody> <tr> <td>FMO 1</td> <td style="background-color: #d9ead3;">-</td> <td style="background-color: #f4cccc;"></td> <td style="background-color: #f4cccc;"></td> <td style="background-color: #f4cccc;"></td> </tr> <tr> <td>FMO 2</td> <td style="background-color: #d9ead3;"></td> <td style="background-color: #f4cccc;">-</td> <td style="background-color: #f4cccc;"></td> <td style="background-color: #f4cccc;"></td> </tr> <tr> <td>FMO 3</td> <td style="background-color: #d9ead3;"></td> <td style="background-color: #f4cccc;"></td> <td style="background-color: #f4cccc;">-</td> <td style="background-color: #f4cccc;"></td> </tr> <tr> <td>FMO 4</td> <td style="background-color: #d9ead3;"></td> <td style="background-color: #f4cccc;"></td> <td style="background-color: #f4cccc;"></td> <td style="background-color: #f4cccc;">-</td> </tr> </tbody> </table>		FITC	PE	APC	Cy7	FMO 1	-				FMO 2		-			FMO 3			-		FMO 4				-	A series of controls with all but one of your colours over every fluorophore (see table) Mandatory for every colour the first time you run your assay. You may not need to run a FMO for every colour on day-to-day runs :D
	FITC	PE	APC	Cy7																							
FMO 1	-																										
FMO 2		-																									
FMO 3			-																								
FMO 4				-																							

Controls for

De-convolutes negative and positive populations within a complex colour assay better than single-stain or isotype controls.
See Detmold staff, or flow-literature for a through explanation.